

---

## Product Data Sheet

---

Product Name: ABM-14  
Cat. No.: GC62391

### Chemical Properties

Cas. No. 1973408-76-2

Formula  $C_{25}H_{18}F_3N_3O_2S$  M.Wt 481.49

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

ABM-14 is a ligand for targeting androgen receptor (AR) for PROTAC. ABM-14 binds to a ligand for VHL via linker to form ARCC-4 to degrade AR[1].

[1]. Meizhong Jin, et al. Compounds and methods for the targeted degradation of the androgen receptor. WO2016118666A1.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---