
Product Data Sheet

Product Name: Ginger extract

Cat. No.: GC62540

Chemical Properties

Cas. No. 84696-15-1

Formula

M.Wt

Solubility DMSO : 100 mg/mL (Need ultrasonic)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Ginger Extract is extracted from Zingiber officinale Roscoe, which has antioxidant potential and regulates STAT6 and TNF- α pathways.

[1] Jedli O, et al. 3 Biotech. 2022 Sep;12(9):191.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA