

---

**Product Data Sheet**

---

Product Name: Bis-PEG3-acid

Cat. No.: GC62672

**Chemical Properties**

Cas. No. 96517-92-9

Formula  $C_{10}H_{18}O_7$  M.Wt 250.25Solubility Storage Store at  $-20^{\circ}C$ 

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Bis-PEG3-acid is a PEG-based PROTAC linker can be used in the synthesis of PROTACs.

[1] Zhao Q, et al. J Med Chem. 2019 Oct 24;62(20):9281-9298.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---