
Product Data Sheet

Product Name: 3-Methylhistamine dihydrochloride

Cat. No.: GC62795

Chemical Properties

Cas. No. 36475-47-5

Formula $C_6H_{13}Cl_2N_3$

M.Wt 198.09

Solubility DMSO : 12.5 mg/mL (63.10 mM; ultrasonic and warming and adjust pH to 10 with NaOH and heat to 60°C) Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

3-Methylhistamine dihydrochloride is an endogenous metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA