
Product Data Sheet

Product Name: FAPI-34
Cat. No.: GC62969

Chemical Properties

Cas. No. 2374782-07-5

Formula $C_{50}H_{57}F_2N_{13}O_{18}$ M.Wt 1166.06

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

FAPI-34 is a fibroblast-activating protein (FAP) inhibitor with favorable pharmacokinetic and biochemical properties. (patent WO2019154886A1).

[1]. Uwe Haberkorn, et al. Fap inhibitor. WO2019154886A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
