
Product Data Sheet

Product Name: IL-17 modulator 4

Cat. No.: GC63019

Chemical Properties

Cas. No. 2446803-65-0

Formula $C_{27}H_{34}N_6O_2$ M.Wt 474.6

Solubility DMSO : 300 mg/mL (632.11 mM; Need ultrasonic) Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

IL-17 modulator 4 is a prodrug of IL-17 modulator 1 . IL-17 modulator 1 is an orally active, highly efficacious IL-17 modulator[1].

[1]. Kevin Neil Dack, et al. Amino-acid anilides as small molecule modulators of il-17. WO2020127685A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA