
Product Data Sheet

Product Name: N-Acetyl-L-arginine dihydrate

Cat. No.: GC63092

Chemical Properties

Cas. No. 210545-23-6

Formula $C_8H_{20}N_4O_5$ M.Wt 252.27

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

N-Acetyl-L-arginine (Ac-Arg-OH) dihydrate is an endogenous metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA