
Product Data Sheet

Product Name: Potassium 2-hydroxy-2-methylsuccinate

Cat. No.: GC63152

Chemical Properties

Cas. No. 1030365-02-6

Formula $C_5H_6K_2O_5$

M.Wt 224.29

Solubility H_2O : 250 mg/mL (1114.63 mM; Need ultrasonic)

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Potassium 2-hydroxy-2-methylsuccinate is an endogenous metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA