

---

## Product Data Sheet

---

Product Name: Biotin-PEG5-azide

Cat. No.: GC63307

### Chemical Properties

Cas. No. 1163732-89-5

Formula  $C_{22}H_{40}N_6O_7S$  M.Wt 532.65

Solubility Storage Store at  $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Biotin-PEG5-azide is a PEG-based PROTAC linker can be used in the synthesis of PROTAC.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA