

---

**Product Data Sheet**

---

Product Name: Mal-PEG4-VC-PAB-DMEA

Cat. No.: GC63400

**Chemical Properties**

Cas. No. 1569261-93-3

Formula  $C_{41}H_{65}N_9O_{13}$  M.Wt 892.01

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Mal-PEG4-VC-PAB-DMEA is a cleavable ADC linker containing a Maleimide group. Mal-PEG4-VC-PAB-DMEA is used in the synthesis of antibody-drug conjugates (ADCs)[1].

[1]. Italo Beria, et al. New alkylating agents. WO2013149948A1.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---