

---

**Product Data Sheet**

---

Product Name: Bersacapavir

Cat. No.: GC63421

**Chemical Properties**

Cas. No. 1638266-40-6

Formula  $C_{16}H_{14}F_4N_4O_3S$ 

M.Wt 418.37

Solubility DMSO : 100 mg/mL (239.02 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Bersacapavir is a novel Hepatitis B Virus capsid assembly modulator[1].

[1]. Vandebossche J, Jessner W, van den Boer M, Biewenga J, Berke JM, Talloen W, De Zwart L, Snoeys J, Yogaratnam J. Pharmacokinetics, Safety and Tolerability of JNJ-56136379, a Novel Hepatitis B Virus Capsid Assembly Modulator, in Healthy Subjects. *Adv Ther.* 2019 Sep;36(9):2450-2462.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA