
Product Data Sheet

Product Name: Atopaxar

Cat. No.: GC63658

Chemical Properties

Cas. No. 751475-53-3

Formula $C_{29}H_{38}FN_3O_5$

M.Wt 527.63

Solubility DMSO : 250 mg/mL (473.82 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Atopaxar(E5555,ER-172594-00) is a potent, orally active, selective and reversible thrombin receptor protease-activated receptor-1 (PAR-1) antagonist.

[1] Dockendorff C, et al. ACS Med Chem Lett. 2012 Mar 8;3(3):232-237.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
