
Product Data Sheet

Product Name: MLKL-IN-2

Cat. No.: GC63680

Chemical Properties

Cas. No. 899759-16-1

Formula $C_{26}H_{25}N_5O$ M.Wt 423.51Solubility Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

MLKL-IN-2 is a MLKL inhibitor extracted from patent WO2021224505A1, compound (i)[1].

[1]. Sáez AJG, et, al. Modulation of mixed lineage kinase domain-like protein signaling. WO2021224505A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
