
Product Data Sheet

Product Name: Propargyl-PEG2-CH₂COOH

Cat. No.: GC63696

Chemical Properties

Cas. No. 944561-46-0

Formula C₉H₁₄O₅ M.Wt 202.2

Solubility Storage 4°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Propargyl-PEG2-CH₂COOH is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs[1].

[1]. Guangrong Zheng, et al. Compounds that induce degradation of anti-apoptotic bcl-2 family proteins and the uses thereof. WO2017184995A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
