

---

## Product Data Sheet

---

Product Name: TLR4-IN-C34-C2-amide-C6-OH

Cat. No.: GC63771

### Chemical Properties

Cas. No.

Formula  $C_{25}H_{42}N_2O_{11}$  M.Wt 546.61

Solubility Storage Store at  $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

TLR4-IN-C34-C2-amide-C6-OH is a linker that incorporates TLR4 inhibitor TLR4-IN-C34. TLR4-IN-C34 inhibits TLR4 in enterocytes and macrophages, and reduces systemic inflammation in mouse models of endotoxemia and necrotizing enterocolitis[1].

[1]. Neal MD, et al. Discovery and validation of a new class of small molecule Toll-like receptor 4 (TLR4) inhibitors. PLoS One. 2013;8(6):e65779. Published 2013 Jun 12.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA