
Product Data Sheet

Product Name: BIIB091
Cat. No.: GC63846

Chemical Properties

Cas. No. 2247614-80-6

Formula C₂₈H₃₄N₁₀O₂ M.Wt 542.64

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

BIIB091 is a highly selective, reversible BTK inhibitor for treating autoimmune diseases.

[1]. Hopkins BT, et al. Discovery and Preclinical Characterization of BIIB091, a Reversible, Selective BTK Inhibitor for the Treatment of Multiple Sclerosis. J Med Chem. 2021 Nov 4.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA