
Product Data Sheet

Product Name: Abiraterone decanoate

Cat. No.: GC64010

Chemical Properties

Cas. No. 2486052-18-8

Formula C₃₄H₄₉NO₂ M.Wt 503.76

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Abiraterone decanoate is a potent Abiraterone prodrug. Abiraterone decanoate provide a controlled release of Abiraterone and long-acting CYP17 inhibition with intramuscular (IM) delivery[1].

[1]. William R. Moore, et al. Abiraterone decanoate (AD): Potent and long-acting activity of a novel intramuscular (IM) abiraterone prodrug depot in a castrate monkey model. Journal of Clinical Oncology, Volume 39, Issue 6_suppl.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA