

---

## Product Data Sheet

---

Product Name: Bis-sulfone-PEG3-Azide

Cat. No.: GC64023

### Chemical Properties

Cas. No. 1802908-01-5

Formula C33H40N4O9S2

M.Wt

700.82

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Bis-sulfone-PEG3-Azide is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs).

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA