
Product Data Sheet

Product Name: MP07-66
Cat. No.: GC64031

Chemical Properties

Cas. No. 1938056-90-6

Formula C₁₉H₃₃NO₃

M.Wt 323.47

Solubility DMSO : 250 mg/mL (772.87 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

MP07-66, a FTY720 analogue, is devoid of immunosuppressive effects and shows promising antitumor effects in chronic lymphocytic leukemia by disruption of the SET-PP2A complex leading to PP2A reactivation.

[1] Zonta F, et al. Blood. 2015 Jun 11;125(24):3747-55.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA