
Product Data Sheet

Product Name: 2-Hexyldecanoic acid

Cat. No.: GC64054

Chemical Properties

Cas. No. 25354-97-6

Formula C₁₆H₃₂O₂

M.Wt 256.42

Solubility DMSO : 100 mg/mL (389.99 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

2-hexyldecanoic acid can be used as an amidating agent [1].

[1]. •Lydie Antonietti, et al. Core-Shell-Structured Highly Branched Poly(ethylenimine amide)s: Synthesis and Structure.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
