
Product Data Sheet

Product Name: AZ14145845

Cat. No.: GC64071

Chemical Properties

Cas. No.

Formula C32H35N9O

M.Wt

561.68

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

AZ14145845 is a highly selective type I1/2 dual Mer/Axl kinase inhibitor with in vivo efficacy.

[1]. McCoull W, et al. Optimization of an Imidazo[1,2-a]pyridine Series to Afford Highly Selective Type I1/2 Dual Mer/Axl Kinase Inhibitors with In Vivo Efficacy. J Med Chem. 2021 Sep 23;64(18):13524-13539.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
