
Product Data Sheet

Product Name: Selnoflast

Cat. No.: GC64089

Chemical Properties

Cas. No. 2260969-36-4

Formula C₂₀H₂₉N₃O₃S

M.Wt 391.53

Solubility DMSO : 12.5 mg/mL (31.93 mM; ultrasonic and warming and heat to 60°C)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Selnoflast (example 6) is a NLRP3 inhibitor (extracted from patent WO2019008025)[1].

[1]. Cooper, Matthew, et al. Novel sulfonamidocarboxamide compounds as NLRP3 inhibitors and their preparation. Patent. WO2019008025.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA