
Product Data Sheet

Product Name: Amino-PEG1-C2-acid

Cat. No.: GC64144

Chemical Properties

Cas. No. 144942-89-2

Formula C₅H₁₁NO₃ M.Wt 133.15

Solubility Storage 4°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Amino-PEG1-C2-acid is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA