
Product Data Sheet

Product Name: Azido-PEG1-azide

Cat. No.: GC64170

Chemical Properties

Cas. No. 24345-74-2

Formula C₄H₈N₆O

M.Wt 156.15

Solubility DMSO : 100 mg/mL (640.41 mM; Need ultrasonic)

Storage -20°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Azido-PEG1-azide is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs[1].

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

[1]. An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA