
Product Data Sheet

Product Name: GPR84 antagonist 1

Cat. No.: GC64229

Chemical Properties

Cas. No. 2762180-28-7

Formula C₂₆H₂₂N₄O₂ M.Wt 422.48

Solubility DMSO : 50 mg/mL (118.35 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

GPR84 antagonist 1 is a high affinity and highly selective competitive antagonist of human GPR84.

[1]. Jenkins L, et al. Discovery and Characterization of Novel Antagonists of the Proinflammatory Orphan Receptor GPR84. ACS Pharmacol Transl Sci. 2021 Sep 7;4(5):1598-1613.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA