

---

## Product Data Sheet

---

Product Name: IL-17A modulator-2

Cat. No.: GC64269

### Chemical Properties

Cas. No. 2748749-47-3

Formula C<sub>33</sub>H<sub>31</sub>N<sub>5</sub>O<sub>5</sub>

M.Wt 577.63

Solubility DMSO : 100 mg/mL (173.12 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

IL-17A modulator-2 is a IL-17A modulator, extracted from patent WO2021239743+A1, example 27. IL-17A modulator-2 inhibits the biological action of IL-17A with a pIC<sub>50</sub> of 8.3. IL-17A modulator-2 can be used for the research of diseases or disorders associated with modulation of IL-17A activity including diseases with an immune component or autoimmune pathol, cancer and neurodegenerative disorders[1].

[1]. Barrie Phillip MARTIN, et al. Preparation of amino acid amides as IL-17A modulators. WO2021239743 A1.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA