
Product Data Sheet

Product Name: BAY 2476568

Cat. No.: GC64302

Chemical Properties

Cas. No. 2311901-93-4

Formula C₂₄H₂₇FN₄O₄

M.Wt

454.49

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

BAY 2476568 is a potent and selective EGFR inhibitor, with IC₅₀s of < 0.2 nM for wild-type EGFR and several mutations (EGFR ex20insSVD, EGFR ex20insASV, EGFR ex20insNPG)[1].

[1]. Stephan Siegel, et al. 4h-pyrrolo[3,2-c]pyridin-4-one derivatives. WO2019081486A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA