

---

**Product Data Sheet**

---

Product Name: Aganepag isopropyl

Cat. No.: GC64435

**Chemical Properties**

Cas. No. 910562-20-8

Formula C<sub>27</sub>H<sub>37</sub>NO<sub>4</sub>S M.Wt 471.65

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Aganepag isopropyl (AGN-210961) is an EP2 agonist[1][2].

Aganepag isopropyl (AGN-210961) is an EP2 agonist[1][2].

[1]. Impagnatiello F, et al. Prostaglandin analogues and nitric oxide contribution in the treatment of ocular hypertension and glaucoma. Br J Pharmacol. 2019;176(8):1079-1089.

[2]. Matsou A, et al. Investigational drugs targeting prostaglandin receptors for the treatment of glaucoma. Expert Opin Investig Drugs. 2018;27(10):777-785.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA