
Product Data Sheet

Product Name: Biotin-nPEG-amine

Cat. No.: GC64454

Chemical Properties

Cas. No. 604786-74-5

Formula C₁₆H₃₀N₄O₄S M.Wt

Solubility Storage 4°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Biotin-nPEG-amine is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA