
Product Data Sheet

Product Name: RapiFluor-MS

Cat. No.: GC64589

Chemical Properties

Cas. No. 1429047-69-7

Formula C₂₁H₂₅N₅O₅ M.Wt 427.45

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

RapiFluor-MS labeling used for LC-MS/MS analysis of N-glycans[1].

[1]. Shiyue Zhou, et al. Direct comparison of derivatization strategies for LC-MS/MS analysis of N-glycans. Analyst. 2017 Nov 20;142(23):4446-4455.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA