

---

**Product Data Sheet**

---

Product Name: BAY 2666605

Cat. No.: GC64703

**Chemical Properties**

Cas. No. 2275774-60-0

Formula C<sub>17</sub>H<sub>12</sub>F<sub>4</sub>N<sub>2</sub>O<sub>2</sub>

M.Wt

352.28

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

BAY 2666605 is an orally active phosphodiesterase 3A (PDE3A) and phosphodiesterase 3B (PDE3B) inhibitor with IC<sub>50</sub>s of 87 nM and 50 nM, respectively.

[1] Manuel E, et al. WO2019025562A1.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA