
Product Data Sheet

Product Name: MC-Val-Cit-PAB-Amide-TLR7 agonist 4

Cat. No.: GC64710

Chemical Properties

Cas. No. 2413016-49-4

Formula C52H72N12O11

M.Wt

1041.2

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

MC-Val-Cit-PAB-Amide-TLR7 agonist 4 (example 15) is a HER2-TLR7 and HER2-TLR8 immune agonist conjugate[1].

[1]. Valerie Odegard, et al. Compositions for the treatment of disease with immune stimulatory conjugates. WO2020056008 A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
