
Product Data Sheet

Product Name: PRMT5-IN-20

Cat. No.: GC65027

Chemical Properties

Cas. No. 880813-30-9

Formula C₂₁H₂₁N₃ M.Wt 315.41

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

PRMT5-IN-20 is a selective protein arginine methyltransferase 5 (PRMT5) inhibitor with anti-tumor activity[1].

[1]. Robert Baiocchi, et al. Compositions and methods for cancer detection and treatment. WO2011079236.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA