
Product Data Sheet

Product Name: Belimumab

Cat. No.: GC65031

Chemical Properties

Cas. No. 356547-88-1

Formula M.Wt

Solubility Storage Store at 4°C, do not freeze

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Belimumab (LymphoStat B) is a human IgG1 λ monoclonal antibody that inhibits B-cell activating factor (BAFF). Belimumab can be used for systemic lupus erythematosus (SLE) research[1].

[1]. Hannah A Blair, et al. Belimumab: A Review in Systemic Lupus Erythematosus. *Drugs*. 2018 Mar;78(3):355-366.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
