
Product Data Sheet

Product Name: PDE11-IN-1

Cat. No.: GC65039

Chemical Properties

Cas. No. 522652-41-1

Formula C₁₆H₁₀CIN₃O₃S

M.Wt

359.79

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

PDE11-IN-1 is a PDE11 inhibitor and can be used for adrenal insufficiency research[1].

[1]. Charles Stuart Hoffman, et al. Inhibitors of phosphodiesterase 11 (PDE11). US9173884.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA