
Product Data Sheet

Product Name: Alirocumab (anti-PCSK9)

Cat. No.: GC65044

Chemical Properties

Cas. No. 1245916-14-6

Formula M.Wt

Solubility Storage Store at 4°C, Do not freeze

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Alirocumab (anti-PCSK9) works by inhibiting the PCSK9 protein, which binds to the low-density lipoprotein receptor (LDLR) and causes receptor degradation, removing less LDL cholesterol from the circulation.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA