
Product Data Sheet

Product Name: Sintilimab

Cat. No.: GC65290

Chemical Properties

Cas. No. 2072873-06-2

Formula M.Wt

Solubility Storage Store at 4°C, do not freeze

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Sintilimab (anti-PD-1, IBI 308) is a fully human IgG4 monoclonal antibody that binds to PD-1 on the surface of T-cells, blocks the PD-1/PD-Ligand 1 (PD-L1) pathway, and reactivates T-cells to kill cancer cells.

[1] Lin Zhang, et al. Front Oncol . 2020 Nov 26;10:594558.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA