
Product Data Sheet

Product Name: PAB-Val-Lys-Boc

Cat. No.: GC65427

Chemical Properties

Cas. No. 1432969-86-2

Formula C₂₃H₃₈N₄O₅

M.Wt 450.57

Solubility DMSO : 125 mg/mL (277.43 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Protocol

Cleavable

Background

PAB-Val-Lys-Boc is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs) for stimulating an immune response[1].

[1]. Geoffrey Alan Pietersz. Immunomodulatory conjugates. WO2013067597A1 (compound PEG-A)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA