

---

**Product Data Sheet**

---

Product Name: Sodium phenyl ethylamido hyaluronate(10% substitution)

Cat. No.: GC65522

**Chemical Properties**

Cas. No.

Formula

M.Wt

Solubility H<sub>2</sub>O : 100 mg/mL (Need ultrasonic)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Sodium phenyl ethylamido hyaluronate (10% substitution) is a chemical agent. Sodium phenyl ethylamido hyaluronate (10% substitution) has 10% phenethylamine substitution. The molecular weight is 40-60WD.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA