

---

## Product Data Sheet

---

Product Name: Methyl heptadecanoate-d33

Cat. No.: GC65526

### Chemical Properties

Cas. No. 1219804-81-5

Formula  $C_{18}H_{33}D_{33}O_2$  M.Wt 317.68

Solubility DMSO : 66.67 mg/mL (209.87 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Methyl heptadecanoate-d33 is the deuterium labeled Methyl heptadecanoate. Methyl heptadecanoate is a fatty acid methyl ester[1].

Stable heavy isotopes of hydrogen, carbon, and other elements have been incorporated into drug molecules, largely as tracers for quantitation during the drug development process. Deuteration has gained attention because of its potential to affect the pharmacokinetic and metabolic profiles of drugs[1].

[1]. Russak EM, et al. Impact of Deuterium Substitution on the Pharmacokinetics of Pharmaceuticals. *Ann Pharmacother.* 2019;53(2):211-216. [2]. Evandro Pereira, et al. Development and validation of analytical methodology by GC-FID using hexadecyl propanoate as an internal standard to determine the bovine tallow methyl esters content. *J Chromatogr B Analyt Technol Biomed Life Sci.* 2018 Sep 1;1093-1094:134-140.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA