
Product Data Sheet

Product Name: CD38 inhibitor 2

Cat. No.: GC65530

Chemical Properties

Cas. No. 2597933-78-1

Formula C₁₉H₂₄N₆O₃ M.Wt 384.43

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

CD38 inhibitor 2 is a potent CD38 inhibitor (IC₅₀ = 0.01 ~ 0.1 μμ).

[1]. Laurie B. Schenkel, et al. Heterobicyclic amides as inhibitors of cd38. WO2021021986 A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
