
Product Data Sheet

Product Name: Aminoxy-amido-PEG4-propargyl

Cat. No.: GC65693

Chemical Properties

Cas. No. 2253965-03-4

Formula $C_{13}H_{24}N_2O_6$ M.Wt 304.34

Solubility Storage 4°C, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Aminoxy-amido-PEG4-propargyl is a non-cleavable 4 unit PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs)^[1].

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA