

---

## Product Data Sheet

---

Product Name: Boc-NH-PEG11-NH2

Cat. No.: GC65724

### Chemical Properties

Cas. No. 1233234-77-9

Formula C<sub>29</sub>H<sub>60</sub>N<sub>2</sub>O<sub>13</sub> M.Wt 644.79

Solubility Storage 4°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

### Background

Boc-NH-PEG11-NH2 is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs<sup>[1]</sup>.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA