
Product Data Sheet

Product Name: N-(t-Boc-Aminoxy-PEG2)-N-bis(PEG3-propargyl)

Cat. No.: GC65725

Chemical Properties

Cas. No. 2112737-60-5

Formula C₂₉H₅₂N₂O₁₁ M.Wt 604.73

Solubility Storage 4°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

N-(t-Boc-Aminoxy-PEG2)-N-bis(PEG3-propargyl) is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs^[1].

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA