

---

**Product Data Sheet**

---

Product Name: Azide-PEG5-Tos

Cat. No.: GC65726

**Chemical Properties**

Cas. No. 236754-49-7

Formula  $C_{17}H_{27}N_3O_7S$  M.Wt 417.48

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Azide-PEG5-Tos (Azido-PEG5-OTs) is a cleavable 5-unit PEG ADC linker applied into the synthesis of antibody-drug conjugates (ADCs).

[1] Alain Beck, et al. Nat Rev Drug Discov. 2017 May;16(5):315-337.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---