

---

**Product Data Sheet**

---

Product Name: m-PEG16-alcohol

Cat. No.: GC65799

**Chemical Properties**

Cas. No. 133604-58-7

Formula  $C_{33}H_{68}O_{17}$  M.Wt 736.88

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

m-PEG16-alcohol is a PEG-based **PROTAC linker** can be used in the synthesis of PROTACs<sup>[1]</sup>.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---