
Product Data Sheet

Product Name: Amino-PEG8-amine

Cat. No.: GC65837

Chemical Properties

Cas. No. 82209-36-7

Formula $C_{18}H_{40}N_2O_8$ M.Wt 412.52

Solubility Storage 4°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Amino-PEG8-amine is a PEG-based (8 units) PROTAC linker can be used in the synthesis of PROTACs.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA