
Product Data Sheet

Product Name: Tr-PEG4-OH

Cat. No.: GC65839

Chemical Properties

Cas. No. 125274-16-0

Formula $C_{27}H_{32}O_5$ M.Wt 436.54

Solubility Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Tr-PEG4-OH is a non-cleavable 4 unit PEG **ADC** linker used in the synthesis of antibody-drug conjugates (ADCs).

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA