
Product Data Sheet

Product Name: Azido-PEG1-CH₂CO₂H

Cat. No.: GC65857

Chemical Properties

Cas. No. 79598-48-4

Formula C₄H₇N₃O₃

M.Wt

145.12

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Azido-PEG1-CH₂CO₂H is a **PROTAC linker**, which refers to the alkyl/ether composition.

Azido-PEG1-CH₂CO₂H can be used in the synthesis of PROTAC BRD4 Degrader-1^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA