
Product Data Sheet

Product Name: Bis-PEG3-NHS ester

Cat. No.: GC65859

Chemical Properties

Cas. No. 1314378-16-9

Formula $C_{18}H_{24}N_2O_{11}$ M.Wt 444.39

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Bis-PEG3-NHS ester is a noncleavable 3-unit PEG linker for antibody-drug-conjugation (ADC).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA